





Bhat Bio-Scan[®] Viral Transport Medium (VTM 3.0 ML)

With 3 ml Viral Transport Medium in a 10 ml self-standing tube recommended for collection and transport of viruses


VTMs For Viral Flu Checkup



Features of the Kit

- ▶ 3ml Transport Medium in a 10ml Self standing tube
- ▶ Flocked Nylon swabs with break point for collection of samples for Higher sample yield.
- ▶ Compatible for molecular based assays.
- ▶ Designed to maintain the viability and the virulence of the viral sample.
- ▶ Antibiotics to control microbial and fungal contamination.
- ▶ Buffers to control the pH
- ▶ Phenol red as a pH indicator
- ▶ Storage Temperature 15 to 30°C and After Sample Collection Store & Transport at 2-8°C
- ▶ Pack Size : 50S

Throat Swab Sample Collection



Nasal Swab Sample Collection

